





**Orthopaedic saw blades
2010**


Dimensions and descriptions



①

Width

The width of the blade measured at the outer edge of the two extreme teeth.




④

②

Depth

The maximum depth to which the blade is able to penetrate the bone




④

③

Cut thickness

The thickness of the blade measured at the cutting edge.



④

④


Body thickness

The thickness of the blade not including the offset tooth design.

MicroAire is a trademark of MicroAire Surgical Instruments.
DeSoutter is a trademark of DeSoutter Surgical.
3M is a trademark of Linvatec Corporation.
Synthes is a trademark of Synthes Inc.
Stryker is a trademark of Stryker Corporation.
Hall is a trademark of Linvatec Corporation.
Aesculap is a trademark of Aesculap AG& Co. KG.
Medicom is not affiliated with any of the above.

All rights of distribution, also by photocopy, reprint(also in the form of extracts) or storing and recovery in any kind of data processing device are reserved and require our written approval.

Numbering system



Feel the innovation! Feel the performance!

Medicom Premium saw blade


<Unique 'wave' surface>

Medicom Premium cut blade has wave surface which is patent pending for it's outstanding design and performance. This unique design provide

- Compare with other same thickness blade, Medicom Premium blade is lighter than others. It helps to prolong lifetime of power tools.
- Due to the wave surface Medicom premium blade has low surface friction and it increases cutting performance.
- Wave surface facilitates to get removed bone chip during cutting the bones.




Model No. 100




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	Order Code					
									A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	1.27mm	1.27mm	95.5mm	25.4mm	PMS-100A	PSO-100A	PLT-100A	PLM-100A	PAO-100A	PDS-100A	PDO-100A	PAS-100A	PPO-100A	PMO-100A
B	1.40mm	1.40mm	95.5mm	25.4mm	PMS-100B	PSO-100B	PLT-100B	PLM-100B	PAO-100B	PDS-100B	PDO-100B	PAS-100B	PPO-100B	PMO-100B

Model No. 110




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	Order Code					
									A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	1.27mm	1.27mm	95.5mm	23mm	PMS-110A	PSO-110A	PLT-110A	PLM-110A	PAO-110A	PDS-110A	PDO-110A	PAS-110A	PPO-110A	PMO-110A
B	1.40mm	1.40mm	95.5mm	23mm	PMS-110B	PSO-110B	PLT-110B	PLM-110B	PAO-110B	PDS-110B	PDO-110B	PAS-110B	PPO-110B	PMO-110B

Model No. 120




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	Order Code					
									A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	1.27mm	1.27mm	95.5mm	19mm	PMS-120A	PSO-120A	PLT-120A	PLM-120A	PAO-120A	PDS-120A	PDO-120A	PAS-120A	PPO-120A	PMO-120A
B	1.40mm	1.40mm	95.5mm	19mm	PMS-120B	PSO-120B	PLT-120B	PLM-120B	PAO-120B	PDS-120B	PDO-120B	PAS-120B	PPO-120B	PMO-120B

Model No. 200




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvecat / 3M	Linvecat / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	25.4mm	BMS-200A	BSO-200A	BLT-200A	BLM-200A	BAO-200A	BDS-200A	BDO-200A	BAS-200A	BPO-200A	BMO-200A	
B	1.20mm	1.20mm	91.5mm	25.4mm	BMS-200B	BSO-200B	BLT-200B	BLM-200B	BAO-200B	BDS-200B	BDO-200B	BAS-200B	BPO-200B	BMO-200B	
C	1.27mm	1.27mm	91.5mm	25.4mm	BMS-200C	BSO-200C	BLT-200C	BLM-200C	BAO-200C	BDS-200C	BDO-200C	BAS-200C	BPO-200C	BMO-200C	
D	1.40mm	1.40mm	91.5mm	25.4mm	BMS-200D	BSO-200D	BLT-200D	BLM-200D	BAO-200D	BDS-200D	BDO-200D	BAS-200D	BPO-200D	BMO-200D	

Model No. 210



Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvecat / 3M	Linvecat / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	23mm	BMS-210A	BSO-210A	BLT-210A	BLM-210A	BAO-210A	BDS-210A	BDO-210A	BAS-210A	BPO-210A	BMO-210A	
B	1.20mm	1.20mm	91.5mm	23mm	BMS-210B	BSO-210B	BLT-210B	BLM-210B	BAO-210B	BDS-210B	BDO-210B	BAS-210B	BPO-210B	BMO-210B	
C	1.27mm	1.27mm	91.5mm	23mm	BMS-210C	BSO-210C	BLT-210C	BLM-210C	BAO-210C	BDS-210C	BDO-210C	BAS-210C	BPO-210C	BMO-210C	
D	1.40mm	1.40mm	91.5mm	23mm	BMS-210D	BSO-210D	BLT-210D	BLM-210D	BAO-210D	BDS-210D	BDO-210D	BAS-210D	BPO-210D	BMO-210D	


Model No. 220



Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvecat / 3M	Linvecat / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	21mm	BMS-220A	BSO-220A	BLT-220A	BLM-220A	BAO-220A	BDS-220A	BDO-220A	BAS-220A	BPO-220A	BMO-220A	
B	1.20mm	1.20mm	91.5mm	21mm	BMS-220B	BSO-220B	BLT-220B	BLM-220B	BAO-220B	BDS-220B	BDO-220B	BAS-220B	BPO-220B	BMO-220B	
C	1.27mm	1.27mm	91.5mm	21mm	BMS-220C	BSO-220C	BLT-220C	BLM-220C	BAO-220C	BDS-220C	BDO-220C	BAS-220C	BPO-220C	BMO-220C	
D	1.40mm	1.40mm	91.5mm	21mm	BMS-220D	BSO-220D	BLT-220D	BLM-220D	BAO-220D	BDS-220D	BDO-220D	BAS-220D	BPO-220D	BMO-220D	

Model No. 230

Model No. 230



Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
															
A	0.97mm	0.97mm	105mm	19mm	BMS-230A	BSO-230A	BLT-230A	BLM-230A	BAO-230A	BDS-230A	BDO-230A	BAS-230A	BPO-230A	BMO-230A	
B	1.20mm	1.20mm	105mm	19mm	BMS-230B	BSO-230B	BLT-230B	BLM-230B	BAO-230B	BDS-230B	BDO-230B	BAS-230B	BPO-230B	BMO-230B	
C	1.27mm	1.27mm	105mm	19mm	BMS-230C	BSO-230C	BLT-230C	BLM-230C	BAO-230C	BDS-230C	BDO-230C	BAS-230C	BPO-230C	BMO-230C	
D	1.40mm	1.40mm	105mm	19mm	BMS-230D	BSO-230D	BLT-230D	BLM-230D	BAO-230D	BDS-230D	BDO-230D	BAS-230D	BPO-230D	BMO-230D	

Model No. 240

Model No. 240

Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	19mm	BMS-240A	BSO-240A	BLT-240A	BLM-240A	BAO-240A	BDS-240A	BDO-240A	BAS-240A	BPO-240A	BMO-240A	
B	1.20mm	1.20mm	91.5mm	19mm	BMS-240B	BSO-240B	BLT-240B	BLM-240B	BAO-240B	BDS-240B	BDO-240B	BAS-240B	BPO-240B	BMO-240B	
C	1.27mm	1.27mm	91.5mm	19mm	BMS-240C	BSO-240C	BLT-240C	BLM-240C	BAO-240C	BDS-240C	BDO-240C	BAS-240C	BPO-240C	BMO-240C	
D	1.40mm	1.40mm	91.5mm	19mm	BMS-240D	BSO-240D	BLT-240D	BLM-240D	BAO-240D	BDS-240D	BDO-240D	BAS-240D	BPO-240D	BMO-240D	



Model No. 250

Model No. 250

Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvac / 3M	Linvac / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	19mm	BMS-250A	BSO-250A	BLT-250A	BLM-250A	BAO-250A	BDS-250A	BDO-250A	BAS-250A	BPO-250A	BMO-250A	
B	1.20mm	1.20mm	91.5mm	19mm	BMS-250B	BSO-250B	BLT-250B	BLM-250B	BAO-250B	BDS-250B	BDO-250B	BAS-250B	BPO-250B	BMO-250B	
C	1.27mm	1.27mm	91.5mm	19mm	BMS-250C	BSO-250C	BLT-250C	BLM-250C	BAO-250C	BDS-250C	BDO-250C	BAS-250C	BPO-250C	BMO-250C	
D	1.40mm	1.40mm	91.5mm	19mm	BMS-250D	BSO-250D	BLT-250D	BLM-250D	BAO-250D	BDS-250D	BDO-250D	BAS-250D	BPO-250D	BMO-250D	





Model No. 260





Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	91.5mm	17mm	BMS-260A	BSO-260A	BLT-260A	BLM-260A	BAO-260A	BDS-260A	BDO-260A	BAS-260A	BPO-260A	BMO-260A	
B	1.20mm	1.20mm	91.5mm	17mm	BMS-260B	BSO-260B	BLT-260B	BLM-260B	BAO-260B	BDS-260B	BDO-260B	BAS-260B	BPO-260B	BMO-260B	
C	1.27mm	1.27mm	91.5mm	17mm	BMS-260C	BSO-260C	BLT-260C	BLM-260C	BAO-260C	BDS-260C	BDO-260C	BAS-260C	BPO-260C	BMO-260C	
D	1.40mm	1.40mm	91.5mm	17mm	BMS-260D	BSO-260D	BLT-260D	BLM-260D	BAO-260D	BDS-260D	BDO-260D	BAS-260D	BPO-260D	BMO-260D	

Model No. 270





Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	85mm	13mm	BMS-270A	BSO-270A	BLT-270A	BLM-270A	BAO-270A	BDS-270A	BDO-270A	BAS-270A	BPO-270A	BMO-270A	
B	1.20mm	1.20mm	85mm	13mm	BMS-270B	BSO-270B	BLT-270B	BLM-270B	BAO-270B	BDS-270B	BDO-270B	BAS-270B	BPO-270B	BMO-270B	
C	1.27mm	1.27mm	85mm	13mm	BMS-270C	BSO-270C	BLT-270C	BLM-270C	BAO-270C	BDS-270C	BDO-270C	BAS-270C	BPO-270C	BMO-270C	
D	1.40mm	1.40mm	85mm	13mm	BMS-270D	BSO-270D	BLT-270D	BLM-270D	BAO-270D	BDS-270D	BDO-270D	BAS-270D	BPO-270D	BMO-270D	

Model No. 280




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvec / 3M	Linvec / Microaire	A/O Synthes Soderm	Order Code	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.97mm	0.97mm	64mm	13mm	BMS-280A	BSO-280A	BLT-280A	BLM-280A	BAO-280A	BDS-280A	BDO-280A	BAS-280A	BPO-280A	BMO-280A	
B	1.20mm	1.20mm	64mm	13mm	BMS-280B	BSO-280B	BLT-280B	BLM-280B	BAO-280B	BDS-280B	BDO-280B	BAS-280B	BPO-280B	BMO-280B	
C	1.27mm	1.27mm	64mm	13mm	BMS-280C	BSO-280C	BLT-280C	BLM-280C	BAO-280C	BDS-280C	BDO-280C	BAS-280C	BPO-280C	BMO-280C	
D	1.40mm	1.40mm	64mm	13mm	BMS-280D	BSO-280D	BLT-280D	BLM-280D	BAO-280D	BDS-280D	BDO-280D	BAS-280D	BPO-280D	BMO-280D	

Model No. 400




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	Order Code					
A	0.7mm	1.05mm	61mm	33mm	BMS-400A	BSO-400A	BLT-400A	BLM-400A	BAO-400A	BDS-400A	BDO-400A	BAS-400A	BPO-400A	BMO-400A

Model No. 410




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	Order Code					
A	0.7mm	1.05mm	85mm	21mm	BMS-410A	BSO-410A	BLT-410A	BLM-410A	BAO-410A	BDS-410A	BDO-410A	BAS-260A	BPO-410A	BMO-410A

Model No. 420




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	Order Code					
A	0.7mm	1.05mm	80mm	17.5mm	BMS-420A	BSO-420A	BLT-420A	BLM-420A	BAO-420A	BDS-420A	BDO-420A	BAS-420A	BPO-420A	BMO-420A

Model No. 430




Code	Body Thickness	Cut Thickness	Depth	Width	Medicom/ Stryker	Stryker old	Linvatec / 3M	Linvatec / Microaire	Order Code					
A	0.7mm	0.96mm	61.1mm	12mm	BMS-430A	BSO-430A	BLT-430A	BLM-430A	BAO-430A	BDS-430A	BDO-430A	BAS-430A	BPO-430A	BMO-430A

Model No. 440





Code	Body Thickness	Cut Thickness	Depth	Width	Order Code									
					Medicom/ Stryker	Stryker old	Lintratec / 3M	Lintratec / Microaire	A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.7mm	0.96mm	30mm	12mm	BMS-440A	BSO-440A	BLT-440A	BLM-440A	BAO-440A	BDS-440A	BDO-440A	BAS-440A	BPO-440A	BMO-440A

Model No. 450




Code	Body Thickness	Cut Thickness	Depth	Width	Order Code									
					Medicom/ Stryker	Stryker old	Lintratec / 3M	Lintratec / Microaire	A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.7mm	0.96mm	31mm	9.2mm	BMS-450A	BSO-450A	BLT-450A	BLM-450A	BAO-450A	BDS-450A	BDO-450A	BAS-450A	BPO-450A	BMO-450A

Model No. 460





Code	Body Thickness	Cut Thickness	Depth	Width	Order Code									
					Medicom/ Stryker	Stryker old	Lintratec / 3M	Lintratec / Microaire	A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.7mm	0.96mm	31mm	6.5mm	BMS-460A	BSO-460A	BLT-460A	BLM-460A	BAO-460A	BDS-460A	BDO-460A	BAS-460A	BPO-460A	BMO-460A

Model No. 470





Code	Body Thickness	Cut Thickness	Depth	Width	Order Code									
					Medicom/ Stryker	Stryker old	Lintratec / 3M	Lintratec / Microaire	A/O Synthes Sodem	De Soutter	De Soutter old	Aesculap	Aesculap old	Medicom old
A	0.7mm	0.96mm	25mm	6mm	BMS-470A	BSO-470A	BLT-470A	BLM-470A	BAO-470A	BDS-470A	BDO-470A	BAS-470A	BPO-470A	BMO-470A

Model No. 700




Code	Body Thickness	Cut Thickness	Depth	Width	Lintratec Hall/ Micro100 Micro Choice	Lintratec Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Lintratec 3M Microdriver	Lintratec 3M Mini driver I	Lintratec 3M Mini driver II
A	0.4mm	0.7mm	34.76mm	16mm							

Model No. 710




Code	Body Thickness	Cut Thickness	Depth	Width	Lintratec Hall/ Micro100 Micro Choice	Lintratec Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Lintratec 3M Microdriver	Lintratec 3M Mini driver I	Lintratec 3M Mini driver II
A	0.4mm	0.7mm	28.11mm	10mm							

Model No. 720




Code	Body Thickness	Cut Thickness	Depth	Width	Lintratec Hall/ Micro100 Micro Choice	Lintratec Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Lintratec 3M Microdriver	Lintratec 3M Mini driver I	Lintratec 3M Mini driver II
A	0.4mm	0.7mm	44.1mm	9.5mm							

Model No. 730




Code	Body Thickness	Cut Thickness	Depth	Width	Lintratec Hall/ Micro100 Micro Choice	Lintratec Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Lintratec 3M Microdriver	Lintratec 3M Mini driver I	Lintratec 3M Mini driver II
A	0.4mm	0.7mm	28mm	8.4mm							

Model No. 740




Code	Body Thickness	Cut Thickness	Depth	Width	Linvac Hall/ Micro100 Micro Choice	Linvac Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Linvac 3M Microdriver	Linvac 3M Mini driver I	Linvac 3M Mini driver II
A	0.4mm	0.7mm	44.1mm	5.5mm							
					SBL-740A	SBA-740A	SBM-740A	SBS-740A	SBV-740A	SBO-740A	SBC-740A

Model No. 750



Code	Body Thickness	Cut Thickness	Depth	Width	Linvac Hall/ Micro100 Micro Choice	Linvac Hall Micro 100 Oscillating	Microaire 7600	Order Code Stryker TPS	Linvac 3M Microdriver	Linvac 3M Mini driver I	Linvac 3M Mini driver II
A	0.4mm	0.7mm	28.11mm	5.5mm							
					SBL-750A	SBA-750A	SBM-750A	SBS-750A	SBV-750A	SBO-750A	SBC-750A

Model No. 600



Code	Body Thickness	Cut Thickness	Depth	Width	Linvac Hall Sturnum	Order Code Linvac 3M Reciprocating	Stryker EHD System 2000
A	0.7mm	1.02mm	70.6mm	12.75mm			
					RLH-600A	RLT-600A	RST-600A

Model No. 610



Code	Body Thickness	Cut Thickness	Depth	Width	Linvac Hall Sturnum	Order Code Linvac 3M Reciprocating	Stryker EHD System 2000
A	0.7mm	1.02mm	70.2mm	9.5mm			
					RLH-610A	RLT-610A	RLT-610A



IMEDICOM Co., Ltd.
#612 Hanlim Human tower, 1-40
Gueungjung-dong, Gunpo-si,
Gyeonggi-do, South Korea 435-050
T. +82 (0)31 479 1156
F. +82 (0)31 479 1157
www.imedicom.co.kr

Presented by